

SALMONELLAE & SHIGELLA

LIM PANG YONG

INTRODUCTION

The presence in foods of any serotype of Salmonella is potentially dangerous as a source of human disease, either directly upon consumption of food, or indirectly through secondary contamination of utensils, processing equipment or processed foods. A further risk arises through induction of the carrier state in food-handlers.

I CULTURE MEDIA*

Nutrient broth

Selenite broth

Tetrathionate broth

Desoxycholate citrate agar (DCA)

Xylose lysine deoxycholate (XLD)

Triple sugar iron agar (TSI)

MacConkey agar (MCA)

GN broth

Salmonella anti-sera : Polyvalent "O" (somatic)

Polyvalent "H" specific and non-specific (flagellar)

* Refer to Appendix B for methods of media preparation.

II APPARATUS

'Waring' blender & flasks

Pipettes

Scissors & forceps

Alcohol (70% v/v) swabs

Plating loops

Inoculating needle

Conical flasks or screw-cap jars, 250 ml

Autoclave

Incubator

Agitated water bath

Weighing balance

Laminar flow chamber

Glass slides

Petri dish (90 x 15 mm)

III SAMPLING PROCEDURE

Refer to "AEROBIC PLATE COUNT" (E-2) Section III.

IV PROCEDURE

A RESUSCITATION (PRE-ENRICHMENT)

1. Weigh 25 g of the above sample and put them into a 'Waring' blender flask and add approximately 225 ml of sterile nutrient broth. Homogenise for 1 min at low speed.

Also blend 50 g of above sample with 200 ml of GN broth for Shigella.

2. Incubate at 35°C for 24 hrs; for Shigella incubate at 35°C for 18 hrs.

B SELECTIVE ENRICHMENT

1. Mix the resuscitated culture gently and add 1 ml each to 10 ml of tetrathionate broth and 10 ml of selenite broth.
2. Incubate the selective enrichment broths at 35°C for 24 hrs.

C PLATING ON SELECTIVE AGAR MEDIA

1. Each culture of enrichment medium is inoculated onto DCA and XLD agar plates. Inoculate the same for Shigella from GN broth culture (from Step A-2) on MCA, DCA and XLD agar plates.
2. Transfer a loopful of culture and streak to obtain isolated colonies.
3. Incubate at 35°C for 24 hrs.
4. Examine the plates for the presence of Salmonella & Shigella colonies.

For Salmonella:

- a) On XLD agar: appear as pink colonies with black centres of H₂S.
- b) on DCA agar: appear as colourless colonies.

For Shigella

- a) On XLD agar: appear as red or pink colour colonies, about 1 mm Ø.
- b) On DCA & MCA: appear as opaque or transparent colonies.

D SCREENING AND BIOCHEMICAL TESTS

1. Pick up a suspected colony with inoculating wire and inoculate the TSI agar slant by streaking the slant and stabbing the butt. Incubate at 35°C for 24 hrs.
2. *Salmonella* cultures typically produce an alkaline (red) slant and acid (yellow) butt, with or without production of H₂S (blackening of butt) in TSI agar. *Shigella* cultures typically produce red slant and yellow butt, with no H₂S or gas.
3. Purify TSI cultures by streaking onto MCA and incubate for 24 hrs at 35°C. Typical colonies appear transparent and colourless, sometimes with a dark centre.
4. Subculture *Salmonella* colony in nutrient broth and incubate at 35°C for 24 hrs. Screen typical *Shigella* cultures in urea agar and motility medium. *Shigella* is urease negative and non-motile.
5. Using the nutrient broth culture as inoculum perform the following biochemical tests:

<u>Tests (Salmonella)</u>	<u>Results</u>	<u>Tests (Shigella)</u>	<u>Results</u>
Lysine decarboxylase	+	Glucose (gas)	-
Urease	-	VP	-
Dulcitol	+	MR	+
KCN	-	Indole	+/-
Malonate	-	Lysine	-
Indole	-	Arginine	+/-
VP	-	Ornithine	+/-
MR	+	Citrate	-
Citrate	+/-	Mannitol	+/-
Lactose	-	Lactose	-
Sucrose	-		

6. Incubate the tests for 24-28 hrs at 35°C.
7. Note that a large percentage of *Salmonella arizonae* strains are negative for dulcitol utilization; positive for malonate and lactose utilization.
8. Perform serological tests for cultures giving reactions typical of *Salmonella* & *Shigella*.


E SEROLOGICAL CONFIRMATION

1. Emulsify the culture in 2 drops of saline on a clean glass slide.
2. Add one loopful of polyvalent "O" antiserum to the first drop only. Use the second drop as a saline control.
3. Tilt the slide back and forth for 1 minute and examine for agglutination. A positive reaction is when there is agglutination in the test mixture but not in the saline control.
4. Repeat similarly with polyvalent "H" antiserum.
5. Salmonella isolates causes agglutination for both antisera.
6. Conduct the serology for Shigella from Step 1 to Step 3.

REFERENCE

Official Methods of Analysis, Association of Official Analytical Chemists, 14th Edition, 1984.

FLOW DIAGRAM OF EXAMINATION PROCEDURES FOR SALMONELLA



FLOW DIAGRAM OF EXAMINATION PROCEDURES FOR SHIGELLA

