

# DETERMINATION OF K-VALUE BY THE FRESHNESS METER

NG CHER SIANG

## INTRODUCTION

The Freshness Meter (Model KV-101, Oriental Co., Tokyo) comprises a water bath incubator chamber fitted with an oxygen electrode, a digital display and a chart recorder. The system uses a series of enzyme reaction steps to break down the ATP and related compounds which are reflected as oxygen consumed during the reactions. The data are presented both in digital read-out and graphically. The ratio of the decrease in oxygen during the various reactions are used to calculate the K-value of the test sample.

The main principles involved in the use of this equipment are presented in Figures 1 and 2.

## REAGENTS

### A: Sample extraction reagents

- a) 10% trichloroacetic acid (TCA)
- b) 10 N potassium hydroxide (KOH)
- c) Methyl red indicator

### B: Reaction reagents

The following reagents are purchased from the equipment supplier or can be mixed in the laboratory.

- a) Reagent B  
Buffer solution
- b) Reagent P is a mixture of alkaline phosphatase and adenosine deaminase.
- c) Reagent E<sub>0</sub> is a mixture of nucleotide phosphorylase and xanthine oxidase.

## PROCEDURE

### A Sample Preparation

Take approximately 2 g fish meat. Add 5 ml of 10% TCA and homogenise. Filter, using a coarse filter paper. Add one to two drops of methyl red indicator, and add 10N KOH to neutralise. This neutralised extract is labelled sample S<sub>1</sub>.

## **B Initial Preparation**

- 1) Set the water incubation bath to 37°C and allow the temperature of the reaction chamber to reach 37°C.
- 2) Place sufficient Reagent B into a bottle and place in water incubation bath. Put in airstone and apply air bubbles to saturate Reagent B with oxygen. Allow Reagent B to equilibrate at 37°C.
- 3) Zero the digital readout and graph recorder by switching the "FUNC" switch to "ZERO". Use the "ZERO" switch to adjust display to  $\pm 0.1$ . Similarly, use the "POSITION" switch on the chart recorder to adjust the pen position to zero on the chart paper. Return the "FUNC" switch to "MEAS".

## **C Preliminary reaction (See Fig. 2)**

- 1) Fill a test tube with 50 ul of Reagent P and 50 ul of extract solution S<sub>1</sub>. Place the test tube in the water bath incubator. Allow reaction to occur for 15 minutes. This reactant solution is labelled S<sub>2</sub>.

## **D Main measurement procedure (See Fig. 2)**

- 1) Add approximately 1.2 ml of preheated Reagent B into the reaction cell.
- 2) Cap the reaction cell tightly. Ensure no air bubbles are trapped inside the reaction cell.
- 3) Switch on chart recorder and stirrer.
- 4) Add accurately 20 ul of sample S<sub>1</sub>. Wait for one minute.
- 5) Add accurately 15 ul of Reagent E<sub>0</sub>.
- 6) Take measurement D<sub>1</sub>.
- 7) Add exactly 40 ul of sample S<sub>2</sub>.
- 8) Take measurement D<sub>2</sub>.
- 9) Computation of results. (See Fig. 1)

## **E Preparing reaction cell for next sample**

- 1) Remove cap and drain out waste reactants.
- 2) Add distilled water into cell and flush several times.

**REFERENCE**

Oriental Electric Co. Ltd. (1987) : Manual of Freshness Tester (Model KV-101)





Fig. 1 Calculation of K-Value

$$K\text{-value (\%)} = \frac{D1}{D2} \times 100$$

**STEP 1**


Original fish extract S<sub>1</sub> is reacted with Reagent, P (alkaline phosphatase + adenosine deaminase)



Resulting solution S<sub>2</sub> contains  $\sum$  HxR + Hx


**STEP 2**

Original fish extract S<sub>1</sub> is reacted with Reagent Eo (nucleotide phosphorylase + xanthine oxidase)



**STEP 3**

Resulting solution from STEP 1, reacted with Eo



**STEP 4**

$$\text{K value (\%)} = \frac{D_1}{D_2} \times 100$$

Fig. 2 Major Steps in using the Freshness Meter