

# COLIFORMS AND ESCHERICHIA COLI

LIM P. Y.

## INTRODUCTION

Coliforms are Gram-negative, non-sporing, facultatively anaerobic rods which ferment lactose, producing acid and gas within 48 hrs and they belong to the family Enterobacteriaceae. The coliform group includes several genera, some of which are of intestinal origin (*Escherichia*) while others are associated with plant and soil material (*Enterobacter*). Thus it is actually a misconception to consider the coliform group as simply an indicator of faecal pollution.

However, generally speaking, it is the count of *E. coli* that is a more reliable indicator of faecal contamination. Its presence indicates recent faecal contamination as it generally does not survive for long in environments other than the intestine.

Faecal coliforms are a group of coliforms capable of fermenting lactose to produce acid and gas at both 37°C and 44.5 ± 0.5°C in 48 hrs and generally contain a high proportion of *E. coli*. As a significant number of non-faecal coliforms can give a positive faecal coliform test, the test can be made more specific for *E. coli* by testing for the production of indole at 44.5 ± 0.5°C.

## I CULTURE MEDIA\*

Brilliant green bile broth (BGB)  
Butterfield's buffered phosphate diluent  
Eosin methylene blue agar (EMB)  
Koser citrate medium  
Lauryl sulphate tryptose broth (LST)  
MRVP medium  
SIM medium  
Simmons citrate agar  
Nutrient broth

\* Refer to Appendix B for methods of media preparation.

## II CHEMICAL REAGENTS#

- a) Kovac's reagent
- b) Methyl red solution
- c) α-naphthol solution (5% w/v)
- d) KOH solution (40% w/v)

# Refer to Appendix D for methods of reagent preparation.

### III APPARATUS

'Waring' blender & flasks	Autoclave
Pipettes	Incubator
Scissors & forceps	Water-bath
Alcohol lamps	Weighing balance
Alcohol (70% v/v) swabs	Laminar flow chamber

### IV SAMPLING PROCEDURE

Refer to 'AEROBIC PLATE COUNT' (E-2) SECTION III

### V SAMPLE PREPARATIONS

Refer to "AEROBIC PLATE COUNT" (E-2) SECTION IV

### VI PROCEDURE

#### A. EXAMINATION FOR PRESUMPTIVE COLIFORMS

1. Select appropriate dilutions and for every dilution, transfer 1 ml aliquots into each of 3 LST tubes.
2. Invert tubes to ensure Durham tubes do not contain gas bubbles.
3. Incubate the tubes at 35°C for 48 hrs.
4. Any tube producing gas is considered positive for the presence of coliforms.

#### B. CONFIRMATION TESTS FOR COLIFORMS

1. Transfer a loopful of suspension from a positive LST tube into a tube of BGB broth.
2. Invert tubes to ensure Durham tubes do not contain gas bubbles.
3. Incubate the BGB tubes at 35°C for 48 hrs.
4. Examine for gas production.
5. Using the MPN Tables (Appendix A), calculate the MPN of coliforms based on the proportion of confirmed LST tubes (with gas production) for 3 consecutive dilutions.

#### C. EXAMINATION FOR PRESUMPTIVE E. COLI

1. Transfer a loopful from each LST tube (with gas production) into a tube of BGB broth, prewarmed to 44.5°C.
2. Incubate the BGB tubes at 44.5°C for 48 hrs.
3. Examine for gas production at 24 hrs and, if negative, again at 48 hrs.
4. Any tube showing gas production is considered positive for the presence of presumptive E. coli.

#### D. CONFIRMATION TESTS FOR E. COLI

1. Subculture all positive BGB tubes by streaking onto plates of EMB agar.
2. Incubate at 35°C for 18-24 hrs.
3. Examine the plates for suspicious E. coli colonies, ie black or dark centred with or without the greenish metallic sheen.
4. Subculture the suspected E. coli colonies in nutrient broth and incubate at 35°C for 18-24 hrs.
5. Perform the following biochemical tests#:  
Indole production  
Methyl-Red & Voges-Proskauer tests  
Citrate utilization  
  
# Refer to Appendix C for biochemical tests procedures.
6. Interpret results as follow:-

	<u>Indole</u>	<u>MR</u>	<u>VP</u>	<u>Citrate</u>
typical <u>E. coli</u>	+	+	-	-
atypical <u>E. coli</u>	-	+	-	-

7. Using the MPN Tables (Appendix A), calculate the MPN of E. coli based on the proportion of BGB tubes in 3 successive dilutions which were shown to contain E. coli.

#### VII CALCULATION OF MPN

$$\text{Most Probable Number (MPN)} = \frac{\text{Index}}{10} \times (450 + W) \times \frac{1}{W}$$

where W : weight of sample in g

Index : from MPN Tables

#### VIII BACTERIOLOGICAL LIMITS OF E. COLI FOR FISH/FISHERY PRODUCTS (COOKED & RAW)

Cooked products : 100 MPN/g

Raw products : 100 MPN/g

#### REFERENCES

- A. Hazzard. (1985). ASEAN Training Course in Fish Quality Control. Training course organised by HAWKAID, Hawkesbury Agricultural College Research and Development Co. Ltd. Chapter: Microbiology in Seafood Quality Control. Section 2:16.  
Chapter: Fish quality control microbiology Section 6:88.

See Reference 2 in E-2.

A flow diagram of the examination procedures for coliforms and E. coli is included as the following figure.

**FLOW DIAGRAM OF EXAMINATION PROCEDURES FOR COLIFORMS AND E. COLI**

